Max Kohl Aktiengesellschaft

Manufacture and Magazine of Physical Apparatus

Chemnitz

Telegraphic address: Physik, Chemnitz. Germany.

➂

ABC - Code 5th Edition used.

Gentlemen,

We herewith beg to recommend ourselves as manufacturers of

Phyiscal Apparatus

and beg you to kindly adress yourselves to our firm, when you contemplate buying such apparatus.

Annexed we have pleasure in sending you some

specimen-pages of our price-list Nº 50

vol -III

which will enable you to judge the great variety of our manufacture.

Our catalogue 50 vol I-III consisting of 1130 pages and containing about 6000 illustrations is sent, free of charge, to all larger schools cultivating education in Physics.

In case no copy of our catalogue should have been sent to you or it should no more be at your disposal, kindly let us know and we shall gladly send you a copy.

Owing to our large stock in physical apparatus we are in a position to execute even important orders in a short delay.


We shall be glad to hear from you soon and to be favoured with your esteemed orders, which will at all times have our best attention.

We are, dear Sirs,

Yours very faithfully

Chemnitz, date of the postmark.
Germany.

Max Kohl Aktiengesellschaft.



63750. 1:16.

63 750. Megadiascope, model A, for projecting diapositives up to 9 × 12 cm and opaque illustrations up to 17 cm diam. projector arc lamp, handregulated, for direct current of 25 to 30 amperes, parabolic mirror 200 mm in diameter cooled with water, achromatic projection objective of 150 mm focal length, focussing being secured by rack and pinion, continuous cooling water trough, objective carrier on slider, diapositive change-frame for diapositives $8^{1}/_{2} \times 8^{1}/_{2}$ cm, $8^{1}/_{2} \times 10$ cm and 9×12 cm, house of oak-wood, projector chamber of metal with air circulation, doors with peep-glasses, stand of pine with case for storing the accessories. The megadiascope is equipped with:

1 Horizontal Projection Apparatus with adjustable plane mirror of first class quality, condenser 150 mm in diameter, projection objective 53 mm in diameter, with tilting erecting mirror silvered on the front in metal mount, and

Megascope for projecting book illustrations, drawings and flat opaque objects by means of reflected light, having first quality illuminating mirror in metal mount, a Voigtländer Heliar of 240 mm focal length and 54 mm aperture, fine focussing being secured by rack

Equipment of the Class Rooms of Public Schools and Colleges.


Lecture Tables for Physical and Chemical Instruction.

We supply the table tops either in **teak** or **oak**. The teak possesses the property of neither cracking, shrinking nor warping in the wet or dry state. Oak which is intended to be worked up into table tops is used only after having lain for a number of years upon the drying heap, and, finally remained in the cut up condition, in our steam drying room for a number of months. Only thoroughly well seasoned wood is employed in our joiners' shops. We utilise pitch pine exclusively for the body of the tables: this wood being quite free from knots and being the most suitable for laboratory purposes. Only the bottoms and inner walls of the drawers are constructed of deal. All tables and cupboards have oak bottom fillets which do not become disfigured when knocked. Unless advised to the contrary, we supply all tops varnished with three coats of linseed oil. If desired, however, the tops can be supplied stained an acid proof black.

All taps, leads, valves, etc. for gas, water, aspirated and compressed air, etc. are constructed in a reliable manner in our own workshops, and only the best of materials are used for the purpose. The gas taps have hose unions bent in an upward direction, thus obviating any tearing or pinching of the hose itself. The leads and piping are also fitted on to the benches in our workshops by workmen possessing years of experience in such work. Our customers are therefore guaranteed the most substantial and reliable workmanship in this connection.

The lecture tables can be made of any length desired; we have built benches up to 18 m in length.

The following lecture tables are those which are intended in the first place to be used in High Schools and Colleges.



50 003. 1:22.


Lecture Table (as suggested by Weinhold), chiefly intended for instruction in Physics, but also adapted for Chemistry in the case of a common room being used for both Physics and Chemistry (W. D., Platte III and Figs 13, 14, 15), see Fig. 50,003.

£ s. d.

	List No. 50,001	50,002	50,003	50,004	50,005
	Length of Table m 3	3.5	4	4.5	5
	Price of Table £21.10.0	24, 4, 0	27. 0. 0	29. 10. 0	32. 0. 0
impregnated with	List No. 50,006 Length of Table m 3 Price of Table £ 24.0.0	50,007 3.5 27. 0. 0	50,008 4 30, 5, 0	50,009 4.5 33. 0. 0	50,010 5 85, 15, 0

The approximate gross weights, with batten packing are: about 330 kg for a length of table of 3 m, 390 kg for 3,5 m, 435 kg for 4 m, 480 kg for 4,5 m and 550 for 5 m.

The table is 90 cm in height and 80 cm in width. Te top is oak, 40 mm thick, is composed of framework and pannellings, and is varnished with three coats of hot linseed oil or stained an acid-proof black (at a slight extra cost). At the left hand side of the table top a slate slab, 54 cm in length and 54 cm in width is let in, upon which work involving the use of acids can be carried out without in



£ s. d.

Prices on Appli-cation

Experimental Switchboards and Experimental Resistances for Educational Institutions and Laboratories.

General.

With the introduction of heavy currents of electricity for teaching purposes, the necessity has arisen for an experimental switchboard in connection with which all experiments can be made which were hitherto carried out in conjunction with galvanic batteries or accumulators. It must be borne in mind that most of the experimental electrical apparatus, especially in regard to the resistance conditions, are constructed in such manner that they can be worked with a few cells, and

with a low voltage require a greater current intensity, while the working voltage is as a rule 220 volts, or

in a number of cases, 110 volts.

On the other hand, the switchboard should permit of such experiments being made which are only rendered possible by heavy currents, such as the demonstration of the electric are light, the singing are, and others.

The switchboard must accordingly be so constructed that on connecting up to a heavy current network, work can be carried on with a low or a high current intensity, as also at a high and a low voltage. The last condition especially must be satisfied in the case of apparatus working with an interrupted current, such as Ampere's table, the spark coil with platinum interrupter, the electric bell, etc. It must be observed that this is not realised simply by inserting resistances in series, since at the moment of breaking the circuit the full network voltage would prevail at the contacts, and the contacts themselves become fused together. Rather must the voltage be distributed by a shunt method of grouping so that it is impossible for the voltage to exceed a certain value.


The switchboard must in addition be provided with measuring instruments for the current and voltage. In the case of the potential, it must be possible to measure the voltage at the apparatus as well as in the network; it is also desirable to measure both the total current taken from the main and the current consumed in the experimental apparatus alone.

TALEDON STATE OF THE STATE OF T

50471. 1:8.

We have constructed a switchboard which entirely fulfils the conditions just prescribed and permits of the carrying out of all experiments which may need to be conducted in teaching or in the laboratory. We have already supplied this type of switchboard to a very large number of educational institutions, and they have met with universal approval, resulting in repart orders being received from a number of towns. The names of the institutions are appended at the end of this section, together with some testimonials appreciative of our switchboards.

A number of forms of switchboards are described in the following pages partly from data obtained from a study of the various voltage conditions and kinds of currents of electricity works and partly as the result of the differing needs and means of the educational institutions.



50735. 1:10.

Prices of Projection Apparatus.


School Projection Apparatus, Model A, Figure 50,735, for projecting apparatus and photographs; with Condenser, Achromatic Objective, Optical Bench, Change Frame with smaller frame inserted for taking diapositives 9×10.5 cm, 8.5×10 cm and 9×12 cm, 1 stage and one movable slider with stand in which to insert the change frame, the stage or other objects; on tall legs.

$ \begin{array}{llllllllllllllllllllllllllllllllllll$	102	122	152
	43	55	65
	150	180	250
With Arc Lamp for hand regulation, No. 50890	50 730	50 735	50 740
	11, 10 0	12, 10, 0	14. 0. 0
With auto-regulating Direct Current Arc Lamp, $No.50898$ $\left\{ egin{array}{lll} List & No. \\ \pounds \end{array} \right.$	50 731	50 736	50 741
	12, 0, 0	13, 0, 0	14. 10. 0
With Direct Current Nernst Projector Namp, No. 50932 $\left\{ egin{array}{lll} { m List\ No.} \\ { m \pounds} \end{array} \right.$	50 732	50 737	50 742
	10, 10, 0	11, 10, 0	13. 0. 0
With Limelight Burner for House gas and Oxygen, No. 50936 $\left\{ egin{array}{ll} { m List\ No.} \\ { m \pounds} \end{array} \right.$	50 733	50 738	50 743
	10. 0, 0	11. 0. 0	12, 10, 0
With Incandescent Spirit Burner, No. 50 964 $\left\{ egin{array}{lll} { t List \ No.} \\ { t \pounds} \end{array} \right.$	50 734	50 739	50 744
	11. 0. 0	12 0. 0	13, 10, 0

Switchboards see page 169, Regulating Resistances for the apparatus with electric arc lamps — see pp. 1226—1228. Fixed Series Resistances, Transformers and flexible triple leads — see pp. 168 and 169. Better Optical and other Outfits — see Nos. 50,750—50,767, p. 161.

The housing is of black iron; it has a door, in the side wall, having dark glass peep holes, and a further observation window on the opposite side. The lamp is manipulated from behind. For this purpose the back is fitted with a sliding door which can be completely removed. The disturbing light radiated from behind is prevented from issuing from the lantern by black curtains, the housing being well ventilated. The condenser is fitted to the front wall, the optical bench being placed in front of the latter, and carrying two sliders: one of those shown in the illustration carries the diapositive holder with change frame, while the other carries the objective-holder together with the objective.

As regards the various systems of lamps, the lantern of this apparatus is arranged so that all the lamps included in the table of prices can be inserted — e. g, instead of the hand-regulated lamp, a lamp with autoregulation, or a Nernst lamp or focus glow lamp, when the candle-power demands are not too great or if a saving in current has to be effected. Gas and spirit incandescent lamps will also fit the lanterns so that they can be



51529. 1:6.

51 533 (51 537), 51 578. 1:8.

Analytical Balance, also suitable for specific gravity determinations, cf. Figure 51,533, with agate planes, gilt or platinised pans, mahogany-glass case with 4 sliding windows, aluminum beam with quick-swing pointer, on white or plack plate glass base, with rider slide when case is closed. The following are not included in the price: the stage, the glass vessel, the Reimann thermometer flask, the rider weights and counterpoises.

List No.	51,531	51,532	51,533	51,534
To carry	5	50	200	1000 grams
Sensitivity	0.1	0.1	0.2	1 mg
Price £	5, 15, 0	6, 10, 0	7, 15, 0	10, 5,0

Prices exclusives of Rider Weights. For Sets of Weights and Rider Weights, see pages 241-243.

- The preceding, without iron or steel parts, agate axes.

List No.	51,535	51,536	51,537	51,538
To carry	5	50	200	1000 grams
Sensitivity	0.1	0.1	0.2	1 mg
Price £	6 5 0	7 0.0	8 10 0	11 0 0


9.10.0


In order to read the tenth-milligrams direct with this quick-swing balance, a power 5 lens-mirror (achromatic objective with silver coating) is fixed on the pillar magnifies, without reflex or distortion, the whole of an opal glass scale, divided in ½,5ths mm, fitted about 20 mm above the ivory scale.

The whole milligrams are read off on the ivory scale; the tenth-milligrams being determined by glancing at the glass scale on the mirror for purposes of confirmation.

As therefore the mirror and the magnification which is read by both eyes is only used at the last moment, the eyes are not in any way strained. The work with both scales is excellent, sure and rapid. The deflection with all loads is, in the mirror, 10 ° and on the ivory scale 2° per milligram.

If Microscope Reading is fitted in place of the mirror the balance is increased in cost by 1. 5. 0






Mechanics of Solids (Statics and Dynamics).

	£ s. d
$51,\!850.$ Collection of Apparatus for demonstrating the Laws of Mechanics, ${ m Figure}$	9. 0.0
The following are included in the above price: 1 Oak Frame, 1 m high and 1 m wide; 1 Inclined Plane; 1 Wedge Apparatus; 2 Weigh Pans; 2 Aluminium Levers with steel axes and steel pins passing through, one firmly fixed, the other arranged for suspending; 1 Set Pulley Blocks with 2 blocks and 3 pulleys arranged in series; 1 Tackle with 3 pulleys; 1 Arbor Wheel with cramp; 2 1 isc-shaped Pendulums of equal length, hanging in various planes, with swivel for adjusting the length; 3 Lenticular Pendulums of 90, 40 and 10 cm length, with swivels for adjusting; 3 Pulleys with cramps (W. D., Fig. 69 A [62 A]; 1 Loose Pulley, of aluminium, with steel axis, bras bow and with 1 hook; 1 Set of Weights with iron weights of 50 g to 5 kg; 1 Set of Double-cheek Weights, each with 20 weights of 100 g and of 50 g.	
51,851. Collection of Apparatus for demonstrating the Laws of Mechanics, Figure, consisting of the following apparatus, these pieces of apparatus being also supplied singly at the prices given: a, b, d1, d2, e, f, g1, h-m	5. 18. 0
(a) 1 Frame, of oack, 1 m high, 1 m wide, with small hooks screwed in (W. D., Fig. 69B [62B])	0. 12. 0
(b) 1 Lever, of aluminium, with steel axis and with steel pins passing through at equal intervals apart, lacquered in two colours, with metal bow, on cramp, Figure	0. 8.0
(c) Pulleys, aluminium, with steel axis and brass bow:	972
1) With 1 hook	0. 2.0
2) With 2 hooks	0. 2.6

Cl. 5713, 5712, 486.





52 241. 1:6. 52 240. 1:6.


	£ s. d.
52,233. Triangular Screw Thread with Nut cut through, Figure, of wood	0. 2.6
52,234. Square Screw Thread, with Nut cut through, Figure, of wood	0. 3.0
52,235. Model of Screw with 4 threads of different shape and pitch, without nuts, of metal, Figure	1. 10. 0
52.236. Simple Machines, Figure. Wedge, Windlass, Capstan, Toothed Gearing with flywheel, Worm Gearing and Screw Jack, mounted on one baseboard	4, 15, 0
Gearing with Screws, Cranks, Links, Rods, Joints and Discs.	
52,237. Fixed Nut with Rotary Screw Spindle, Figure, with angular, accurately turned thread, constructed entirely in iron	1. 0.0
52,238. Fixed Screw Spindle, with Rotary Progressive Nut, Figure	1. 4.0
$52,239$. Screw Spindle, rotary, resting immovably in its bearings, with straight guided, progressive nut, ${ m Figure}$	1. 10. 0
52,240. Rotary Nut, incapable of lateral movement, with straight guided, progressive screw spindle, Figure	1. 10. 0
52,241. Screw with Right-and-Left-Handed Thread, Figure, with nuts	2. 10. 0

Apparatus for showing the formation and mode of action of the screw, see No. 51,905-51,910, p) 273.


Cl. 518, 519, 3786, 521, 3787, 5805, 3788, 520, 3789.








52465. 1:4




52 463. 1:7.




52 466. 1:8.

THE PART HERE IN THE PART OF T	A S. d
52,460. Flat Sheet Iron Box, for storing the rubber cushion (W. I)., page 148)	0. 8.0
52,461. Hydraulic Press , Figure (W. D., Fig. 107 [96]), massively constructed for iron and steel; valves fitted, so as to be detachable, in special valve chambers, giving 3000 kg pressure, with safety valve for 60 atmospheres	12. 0.0
The ram has a cross-section of 50 sq. cm; each atmosphere of pressure thus corresponds to a 50 kg pressure. The press has an efficiency of 60 ats., which is equivalent to 3000 kg pressure.	
52,462. — idem, with Pressure Gauge	14. 10. 0
52,463. Hydraulic Press, of metal, for a pressure of 1500 kg, Figure	8. 10. 0
$52,\!464.$ — $idem,$ with Pressure Gauge	11. 0.0
52,465. Hydraulic Press with Glass Cylinder and visible valves, Figure, with safety valve for 25 Atm., for preventing bursting of the cylinder, for 500 kg pressure	8. 0.0
52,466. — idem, with Pressure Gauge, Figure	

Cl. 707, 709, 708, 5119.



reading to $1/10$ th mm and movable with rack; reading of vernier and level facilitated by mirror illumination. The thermometer is in direct contact with the barometer tube and can be read from the outside. The instrument is suspended on a hook and can be centred	
at the lower end	
52,787. — idem, with tube 12.5 mm wide, Figure	10. '0.
52,788. Travelling and Altitudinal Barometer, after Fortin, Figure, with tube 10 mm width, vernier giving 1/10 th mm, with thermometer, Stand and Universal Suspension, also leather case	7. 10.
52,789. — idem, Station Barometer, on wood board, Figure	




52 907. 1:9.

52 911. 1:10.

52,908. Vacuum Pump with Oil-packed Piston and oil non-return valve, Kohl's system Patent; with Hand Lever and with Plate 240 mm diameter and Baromet		
200 mm high		11. 0.0
52,909 idem, with Fly-wheel, for Hand and Power Drive, without Plate, Figure	е	10. 0.0
52,910. — idem, with Fly-wheel and with Plate 240 mm de diameter and Barome 200 mm high	ter Gauge	12. 10. 0
Largest Receiver for use with above: No. 53,013		
52,911. — idem, driven by a 1/eth H. P. Electric Motor with worm gearing, motor for Direct Current, with Starter, without Plate, Figure		


Vacuum Pumps No. 52,903—52,910 are specially intended for rapid work in laboratories.

Cl. 943, 3652, 4711, 4713.






53 577 (53 578). 1:8.




5**3** 5**83.** 1:11.

	53,576. Vibration Microscope, after Weinhold, Figure (W. D., Fig. 239 [225]), can be used at same time as a Töpler Vibroscope	
**	53,577. 2 Massive Tuning Forks on Stand for the tone $c_0 = 128$ compound vibrations $ut_2 = 256$ v. s.); the pitch of one can be varied by filling with mercury; both with steel mirror and large resonator, Figure (Pogg. Ann. d. Phys. u. Chemie, Vol. 157, 1876, p. 621), for producing any phase differences and impact tones, also for Lissajous's Figures One of the two tuning forks has both limbs bored out and is provided at the yoke with a screw press by means of which mercury can be forced into te limb. In this manner the frequency can be altered by 8 compound vibrations (= 16 v. s.). The resonator of each fork is provided with a gap which can be	27. 10. 0
	increased or decreased by means of a screw so as to suit its tone to that of the fork. The forks have electromagnetic drive.	
	53,578. Variable Tuning Fork $c_0=128$ compound vibrations alone (ut ₂ = 256 v. s.), Figure, on stand and with variable resonator	19. 5.0
	53,579. 2 Tuning Forks with Resonators, same pattern as No. 53,577, but smaller, for the tone $\mathbf{c_1} = 256$ compound vibrations (ut ₃ = 512 v. s.)	21. 0.0
	53,580. Variable Tuning Fork $c_1=256$ compound vibrations (ut ₈ = 512 v. s.), alone, on stand with variable resonator	15. 10. 0
	The state of the s	


Cl. 1172, 5814, 1173, 3381.



53 823. 1:6.



53,822. Large Reflecting Goniometer with circle 230 mm diameter, with concealed graduation | £ s. d. on silver, objective aperture 31 mm, reading by magnifying glass with 2 verniers for The circle and observing telescope have independent motion about the central axis, with screw motion. the alidade along with the telescope, or by clamping the alidade with the telescope and rotating the circle with the crystal. The instrument has two oculars and a number of slots for the collimator; the observing telescope is provided with an extra objective, thus facilitating the focusing of the objects under test. Illustration on application. The measurements can be taken either by firmly clamping the circle carrying the crystal, and by turning 53.823. Reflecting Goniometer, can also be used as a Spectroscope, Figure, with circle 150 mm diameter, concealed, on silver, with reading by means of magnifying glass, to 20 seconds; objective aperture 27 mm. Without Centering Apparatus 27. 10. 0 53,824. Reflecting Goniometer for the Laboratory, with circle 150 mm diameter, graduated in 19. 0.0 53,825. Centering Device for the crystals, fitting the two preceding goniometers. 3. 15. 0 53,826. Reflecting Goniometer, smaller, without magnifying glass reading, Figure, practical model for laboratories; telescope of 19 mm aperture, graduated in 1/2 and with vernier 13. 0.0 53,827. - idem, Figure, with concealed Divided Circle and 2 magnifiers for reading, without 16. 0.0 This pattern can be highly recommended as the graduation cannot be touched by the fingers.



54 382. 1:5.



54 380. 1:5.




54 385. 1:6.




54 386, 54 386 b, 54 386 d, 54 386 f, 54 386 h. 1:5.


£ s. d
4. 0.0
6. 0.0
6. 12. 0
1. 10. 0
1. 6.0
8. 0.0
12. 10. 0









55 262. 1:9.




55 256. 1:9.



55 259. 1:9.




55 261. 1:7.





55 263. 1:7.

55 262. 1:9.		,		25.400	e de							£ s. d.	
55,256. Sectional Model of a Cylinder with simple Slide Valve,	F	iε	ŗ u	re,	0	f v	v 00	d.					
55,256a idem, of iron, smaller, about 3 the size													
55,257 idem, with Farcot Slide Valve, Figure, of wood .													
55,258. — idem, of iron, smaller, about 3/5 the size													
55,259 idem, with Meyer Compound Slide Valve, Figure,													
55,260. — idem, of iron smaller, about 3/5 the size													
55,261 idem, with Rider Slide Valve, Figure, of iron	4	•	*		8 3					÷		12. 0. 0	
55,262 idem, with Woolf Slide Valve, Figure, of wood .													
55,262 a. — idem, of iron, smaller, about 3/5 the size													
55,263. — idem. different pattern, Figure, of iron												THE RESERVE OF THE PROPERTY OF THE PARTY OF	

4899, 1893 3891











60 097. 1:4.

60 098. 1:5.

60 100. 1:5.








60 101. 1:5.

60 103. 2:5.

60,098. Dip Circle, can also be used as galvanoscope, Figure, with carefully balanced magnetic needle, column turning in base, divided circle, movable on hinge	£ s. d. 2. 8.0
60,099. Dip Circle, Figure 60,064, p. 787, can also be used for the fundamental magnetic experiments (see under No. 60,064)	5. 10. 0
60,100. Dip Circle, Figure, with horizontal and vertical circle, on tripod with levelling screws, the needle works in agate bearings. The vertical circle is movable, being 110 mm in diameter	3. 10. 0
60,101. — idem, larger, Figure, vertical circle 150 mm, horizontal cercle 100 mm diameter, with round spirit level	7. 10. 0
60,102. Dip Circle, large type, Figure (M. P., 9th edn., III, Fig. 35; GanMan. Fig. 611; GanRein. Fig. 703), with detachable needle working in carnelian bearings, vertical circle 190 mm diameter. The instrument gives accurate data	12 .10. 0
60,103. Dip Circle, Figure, with vernier reading for the circle and with micrometer adjustment of the vertical circle; diameter of upper circle: 180 mm	11. 0.0



60425-60430. 2:7-1:7.

60407-60416. 1:8-1:27.

Self-Exciting Influence Machine, as suggested by Töpler, with 1 fixed and 1 rotating plate, Figure, with massive ebonite pillars on which the discharging combs are placed, in order to be able to remove the plates forth; with stand of polished mahogany.

List No. 60,397 60,398 60,399 60,460 60,401 60,402 60,403 60,404 60,405 60,406

Diameter of rotating disc. 26 31 36 41 47 52 57 62 75 90 cm


£ 3.10.0 4.10.0 6.0.0 7.5.0 9.5.0 10.10.0 15.15.0 20.0.0 27.0.0 30.0.0

Self-Exciting Influence Machine, as suggested by Töpler, Figure, with 2 fixed and 2 rotating plates, base and pillars of mahogany, discharging combs fixed on ebonite pillars, very neatly constructed and producing twice the quantity of electricity as those with only two plates.

List No. 60,407 60,408 60,409 60,410 60,411 60,412 60,413 60,414 60,415 60,416

Diameter of rotating plates 26 31 36 41 47 52 57 62 75 90 cm
£ 6.15.0 8.10.0 11.0.0 13.10.0 17.5.0 19.0.0 25.0.0 29.5.0 32.10.0 41.5.0


Cl. 5707, 5599, 2184, 5709.



62 299, 62 302. 1:7.




62 306. 1:4.




62 308, 1:7,5.


62,299. Small Half-Ring Electromagnet after du Bois, Figure, with one pair cylindrical and one pair conical pole pieces, one distance piece for I mm pole distance of the conical pole pieces, Stand and Spanner (Ztschr. f. Instrumentenkunde, 1899, p. 363) This magnet is constructed in similar manner to No. 62,292 but is of only half the linear dimensions. The winding is arranged for 8 amps, and when the coils are in series the resistance is approximately 4 ohms; the apparatus can therefore be conveniently connected up to 32 volts.	£ s. d.											
Accessories for preceding Electromagnet.												
62,300. One Pair Cylindrical Pieces for making additional pole pieces, of dynamo steel	1.10.0											
62,301. Universal Slider, cf. Fig. 62,294	4.15.0											
62,302. v. Waltenhofen's Pendulum	3. 0.0											
62,303. Flat Induction Coil, cf. Fig. 62,296	3. 0.0											
62,304. Polarisation Apparatus with divided circle and Faraday's Glass Block	7. 5. 0											
62,305. Holder for a Bismuth Spiral	0.10.0											
62,306. Bismuth Spiral after Lenard, Figure, for determining the Field Density of Magnetic Fields by the change in the Conducting Resistance which bismuth undergoes in the field; in case, with calibration curve; approx. 20 mm diameter												
62,307. — idem, 5 to 6 mm diameter	3. 0.0											
62,308. Electromagnet after Prof. P. Weiss, with coils for water cooling which can be loaded twice as much as in the case of air cooling; Figure; diameter of pole pieces 100 mm; with one pair cylindrical and one pair conical pole pieces; on stout wood board The pole pieces can be adjusted by a micrometer motion. The coils are arranged as near as possible to the pole pieces and the cross section of the iron increases steadily from the pole pieces up to the centre of the magnetic circuit, so that leakage is avoided as much as possible and a very high degree of magnetization attained. Unless otherwise requested the winding is dimensioned for 120 volts. A great advantage lies in the fact that the entire apparatus remains quite cool.												
62,309. — idem, arranged to rotate on iron base so that the magnet can easily be placed in any position, with horizontal divided circle	105. 0. 0											
62,310. Electromagnet after Prof. Weiss, as No. 62,308, but with pole pieces 92 mm diameter	73.10.0											
62,311. — idem, on tripod, as No. 62,309	85. 0.0											




62 383. 1 ; 6.



62 387. 1:5.




62 382. 1:8.



62 386. 1:10.

62,382. Demonstration Apparatus for the Magnetic Field of a horse-shoe Electromagnet, after Götze, Figure	£ s. d. 1. 10. 0
62,383. Ampere's Stand, Figure, reliable pattern, with separate stands; main stand on metal tripod with levelling screws, rotary and fixed in position by screw; wire figures of aluminium with silver points. The following appertain to above: 1 solenoid (W. D., Fig. 536 [509]), 1 simple rectangle (W. D., Fig. 532 [505]), 1 double (astatic) rectangle	=
(W. D., Fig. 534 [507]), 1 wire stirrup with wire coil (W. D., Fig. 535 [508]), 1 boxwood bowl with 2 concentric rings, 1 boxwood bowl with partition wall. Also adapted for experiments as suggested in GanMan., Figs. 739, 740, 746, 751; GanRein., Fig. 891, 893, 877	4. 4. 0
This pattern enables the experiments to succeed with ease. The primary condition is that absolutely pure and dry mercury should be emp oyed. This we supply if desired at a charge of £ 0.8.0. per kilo. A current density of $15-20$ amps. is necessary. Precise instructions for use are given in.	
62.367. Wire Coil Ampere's Stand (W. D., Fig. 537 [510]), Figure on p. 962	0. 8. 0

2673, 3879,




Induction Coil without Interrupter, Type A, B, C, D and E, on Table Stand.


Spark-		For	Wehnelt	Interru	ipter	Fo	r Wehnelf Interr		For Interr	Spark-		
length with Wehnelt- Inter- rupter	Size	Without Primary Pach	Condenser, Coil with ytrope, c. 62533	Type B Without Condenser, Primary Coil with Plug Commutation, cf. Fig. 62513		With C Primary Pach	ondenser, Coil with ytrope, g. 62533	Type D With Condenser, Primary Coil with Plug Commutation, cf. Fig. 62513		With Co Primary Comm cf. Fig	length with the Motor Inter- rupter	
mm		List No.	£ s. d.	List No.	£ s. d.	List No.	£ s. d.	List No.	£ s. d.	List No.	£ s. d.	mm
				1	. Standar	d Secon	dary W in	ding.				11
_	2				-	_	-	-	-	62559	8. 0.0	90
_	3				-		_	_	-	62 560	10. 0.0	120
-	4		-	222	_	-		_	-	62561	12. 0.0	150
250	5	62491	16.10.0	62:508	15,10.0	62 52 5	18.10.0	62:542	17,10,0	62 562	15. 0.0	200
300	6	62492	20, 0.0	62509	19, 0.0	62526	22. 0.0	62 543	21, 0,0	62563	18.10.0	250
350	7 a	62493	26, 0.0	62510	25. 0.0	62527	28.10.0	62544	27.10.0	62564	25. 0.0	300
350	7 b	62494	30, 0.0	62511	26, 0.0	62 528	32.10.0	62 54 5	28,10,0	-	-	300
400	8	62495	37, 0.0	62512	33, 0.0	62 529	39.10.0	62546	35,10,0	62 566	32,10,0	350
450	9	62496	COMMENT - COMMENT	62513	40, 0.0	62 530	47. 0.0	62 547	43, 0.0	62567	40, 0,0	400
500	10	62497	54. 0.0	62514	50, 0.0	62 531	57. 0.0	62548	53. 0.0	62568	50, 0.0	450
550	11	62498	65, 0.0	62515	60, 0.0	62532	68.10.0	62549	63,10,0	62569	60, 0.0	500
600	11a	62499	75. 0.0	62516	70, 0.0	62 533	78.10.0	62550	73,10,0	62570	70. 0.0	550
650	12	62 500	85, 0.0	62517	80, 0.0	62 534	88.10.0	62551	83,10.0	62571	80. 0.0	600
700	12a	62 50 1	95. 0.0	62518	90, 0.0	62535	99. 0.0	62552	94, 0.0	62572	90, 0,0	650
750	13	62 502	105. 0.0	62519	100, 0.0	62536	109. 0.0	62553	104, 0,0	62573	100, 0.0	700
800	13a	62 503	115. 0.0	62 520	110, 0.0	62 537	119. 0.0	62554	114. 0.0	62574	110. 0.0	750
850	14	62504	125. 0.0	62 521	120, 0.0	62538	129, 0.0	62 555	124. 0.0	62575	120, 0.0	800
900	14a	62505	135. 0.0	62 522	130, 0.0	62539	140. 0.0	62 556	135, 0.0	62576	130, 0.0	900
950	15	62506	150. 0.0	62 523	145, 0.0	62540	155. 0.0	62 557	150, 0.0	TO THE PROPERTY OF WAY	145. 0.0	1 10 - W. P. C. W. W. P.
1000	15a	62 507	165. 0.0	62 524	160, 0.0	62 541	170. 0.0	62558	165, 0,0	62578	160. 0.0	1000
				II. L	ow Resis	tance S	econdary	Winding				
350	RJI	62579	34. 0.0	62581	30. 0.0	62583	36.10.0	62585		62587	29,10,0	7//00000
400	RJII	62580	41.10.0	62 582	37.10.0	62584	44. 0.0	62586	40, 0.0	62588	37. 0.0	400
co			a for We			r have	double, 7	b—10,	quadruble	e, and 1	1—15 a,	sextuple £ s. d.
						wall brac	ked, for siz	es 2-5	ur ge gorgene			1. 0. 0
02,009. 3	parking	for size	6_0 *	miniy h	ted on the	wall blac						1, 16, 0
02,090. – eo so :	idem	, for sizes	10 19			8 8 8			\$ \$ \$ WE		1 10/15	2. 10. 0
		, ,, ,,	10-12 .							25 5 6		
2,592	- luem	1 11 11	164-194								10 (0) 700	- 3000 0000

Chemnitz, Germany

Dynamos

for Generating Direct, Alternating and Three-phase Current; also Accessories.





62767. 1:8.


62777. 1:12.

Shunt Dynamos for Direct Current, cf. Fig. 62,774; also constructed as Universal Dynamos, with 4 slip rings, for Direct, Alternating and Three-phase Current, with Drum Armature (cf. Gan.-Man., Fig. 806).

2.70	. D. C. put	Approx.	Single- phase A. C. voltage	Three- phase pres- sure	Ap- prox. power re-		a Steel ire	98(5)	eeds Lamps	Feeds Arc Lamp	Driving	hout Stand, 62774	Driving cf. Figs and 6	62777
Volts	Amps.	арр	approx. pr	prox. Volts	quired HP	Length mm	Diam.	No.	No. Hefner N	No.	List No.	£ s. d.	List No.	£ s. d
					(a) For	Direct	Curr	ent on	ly.				
20	3	3000		_	0,25	120	0,15	4	3	62815	62770	4, 10, 0	62777	6. 0.0
35	4	2800	-	_	0,33	250	0,2	5	8	62816	62771	8, 0,0	62778	10. 0.0
50	4	2200	_	-	0,5	300	0,2	5	12	62816	62772	10. 0.0	62779	12.10.0
65	5	2200	-	827	0,8	400	0,3	7	16	62816	62773	12, 10, 0	62780*	22, 10, 0
(b) For	Direct	Current	, Mon	ophase	A. C.,	Three-	phase	Curre	nt and	Dissymet	rical Tw	phase	A. C.
20	3	3000	14	12	0,25	120	0,15	4	3	62815	62774	5, 15, 0	62781	7. 5.0
35	4	2800	24	20	0,33	250	0,2	5	8	62816	62775	9. 5.0	62782	11. 5.0
50	4	2200	34	30	0,5	300	0,2	5	12	62816	62776	11, 10, 0	62783	14. 0.0
65	5	2200	44	38	0.8	400	0.3	7	18	62816	62776a	14. 0.0	62783a*	24. 0.0


The Machines are of the two-fold type, provided with Drum Armature and are also suitable for continuous

^{*} The Driving Stand for Nos. 62,780 and 62,783 a is constructed in accordance with Fig. 62,783 b. Cl. 6045, 6046.




62 902 A. 1:3.


62 902 B. 1:3.



62 902 C. 1:3.



62902D. 1:3.



62902 E. 1:3.

62,901. Lantern Slides for Explaining the Connections and Mode of Action of Generator No. 62,896, of A. C. Motor Model No. 62,898 and of the Demonstration Apparatus for rotary field experiments, No. 62,902 (W. D., Figs. 579—581, 583—587, 596 and 598) Each

(a) Alternator (Fig. 579),

(b) Magnetic Axis of the Alternator (Fig. 580),

((c) Stator Winding of the Three-phase Machine I (Fig. 581),

d) Curve for Mono-phase A. C. (Fig. 583),
(e) Curve for two-phase A. C. (Fig. 584),
(f) Curve for three-phase Current (Fig. 585),

(g) Star and Delta Connection (Fig. 586),

(h) Stator Winding of three-phase machine II (Fig. 587),
 (i) Connections used as three-phase motor (Fig. 596),

(k) Connections used as mono-phase motor (Fig. 598).

62,902. Demonstration Apparatus for Rotating Field Experiments (Battery A. C. Apparatus), after Weinhold, Figs. A—E (Ztschr. f. d. phys. u. chem. U., 6, p. 7; W. D., Fig. 582).

The complete apparatus consists of: current distributing apparatus, Fig. A; Ferraris Double Coil Fig. B; Ferraris triple coil, Fig. C; 4-coil iron ring, Fig. D; 6-coil iron ring, Fig. E; magnetic needle and sheet-iron disc with agate cap; squirrel-cage rotor; glass disc with iron filings.


The apparatus permits of carrying out the Alternating and Three-phase Current experiments in a convenient and practical manner. The following pertain to the apparatus: 1—2 Bunsen cells or an accumulator the current of which is transformed into Alternating or Three-phase Current by the distributing apparatus, Fig. A. This current flows through double or triple crossed different-coloured wire rings, Fig. B or C, in which a magnetic needle or an iron disc can be inserted. The bodies inserted in the coils rotate under the influence of the rotating field. Instead of these iron rings an iron ring can be attached consisting of four or six coils and wound with covered wires of different colour (Tesla Ring), Figs. D and E; with these it is possible to obtain synchronous rotation by inserting a magnetic needle or an iron disc; or asynchronous rotation may be produced by an iron ring wound with copper wire (short-circuited rotor). It is possible to show the migration of the lines of force of the field in a very pretty manner by means of a glass disc over which iron filings have been strewn.

6. 0.0

0. 1.6

5. 10.

0. 15. 0



Offices and Manufacture.